

All advantages at anderson-negele.com



50019 / 3.3 / 2021-03-26 / SD

Product Information TFP-641, -642, -661, -681, -841, -842, -861, -881

**PHARMA** 

# Temperature Sensor PHARMadapt EPA



#### **Application / Specified Usage**

- · Developed for applications in pharmaceutical industy and biotechnology
- · Temperature measurement especially at small pipe diameters
- · In connection with build-in system PHARMadapt EPA suitable for pipes DN 10...100

## **Authorisations**





#### **Application Examples**

- · Process monitoring
- · Monitoring of CIP-/SIP-cleaning

#### · Hygienic and easy sterilizable installation by using Negele build-in system **PHARMadapt EPA**

- · CIP-/SIP-cleaning up to 140 °C (284 °F)
- · All product contacting materials compliant to FDA
- · Sensor completely made of stainless steel
- Sealing ring according to USP Class VI
- · Compliant to 3-A Sanitary Standard for DIN 11866 series A with DN ≥ 25,

## **Hygienic Design / Process Connection**

DIN 11866 series B with DN ≥ 20, DIN 11866 series C with DN ≥ 1"

#### Temperature sensor TFP-641 with build-in system EPA-18



#### **Features**

- · For small pipe diameters from DN10
- · Easy demounting for cleaning and calibration by clamp system
- · Short reaction time, very compact measure point with leackage control
- · Integrated transmitter available
- · Light weight sensor head, non-sensitive to vibrations
- · Hygienic lid design
- · Electrical connection via M12 plug
- · Material 1.4435 (AISI 316L), material certificate 3.1 in scope of delivery (for all product contacting parts)
- · Quick and easy to install with an orbital welding machine

#### Temperature sensor TFP-661 with transmitter MPU-M and **TAG-number plate**



#### **Options / Accessories**

- · 2x Pt100, optional (not retrofittable) 2x Pt100 with two transmitters (not retrofittable)
- · Programmable transmitter MPU-4 and MPU-M with 4...20 mA output, 2-wire
- · Integrated transmitter for HART-protocol
- · Programming adapter MPU-P 9701
- · Integrated display MPU-LCD in connecting head
- · Pt100-chip with other classes of accuracy (A, AA, AAA)
- · Preassemled cable for M12 plug
- · Fixed cable in other length or material available
- · Calibration certificate (just available with order placement)
- · Customer specific label and TAG-number plate (stainless steel)

Temperature sensor				
Process connection	gap free	with clamp-ring SRC-05 resp. SRC-10		
Insertion length EL	TFP-6xx TFP-8xx	10 mm, 25 mm, 50 mm, 100 mm 20 mm, 50 mm		
Materials	connection head protection tube sealing ring	stainless steel 1.4301 (AISI 304) stainless steel 1.4435 (AISI 316L) EPDM, USP Class VI, FDA 21CFR177.2600		
Temperature ranges	ambient sensor tip	-5080 °C (-58 °F176 °F) -50250 °C (-58 °F482 °F)		
Operating pressure		10 bar max.		
Sensing resistor	acc. to DIN EN 60751	Pt100		
Electrical connection	TFP-64x, TFP-84x  TFP-661, TFP-861  TFP-681, TFP-881	cable gland M16x1.5 (PG) or M12 plug 1.4301 (AISI 304), 4-pin M12 plug 1.4301 (AISI 304) fixed cable (PTFE, 4x 0.14 mm²), standard: 2.5 m		
Protection type		IP 69 K (with electrical connection M12 plug)		

Transmitter MPU-4, MPU-H, MPU-M				
Temperature ranges	ambient storage	-4085 °C (-40185 °F) -5590 °C (-67194 °F)		
Measuring ranges	MPU-4, MPU-H, MPU-M	standard: -1040 °C, 050 / 100 / 150 / 200 °C special ranges free programable		
Accuracy	input	< ±0.25 °C		
Temperature drift	zero, span	< 0.01 % / K		
Supply	MPU-4, MPU-H, MPU-M accuracy	835 V DC 0.01 % / V (reference: 12 V DC)		
Output	signal accuracy burden	analog 420 mA < $\pm 0.1$ % of measurement range < $600 \Omega$ (at U <sub>B</sub> = 24 V)		
Humidity	without condensation	098 %		

Accuracy classes of temperature sensors   Tolerances for Pt100 acc. to DIN EN 60751					
Pt100	Class A	Class AA / 1/3 DIN B	Class AAA / 1/10 DIN B		
0°C/100Ω	±0.15 K / ±0.06 Ω	±0.10 K / ±0.04 Ω	±0.03 K / ±0.01 Ω		
100 °C / 138.5 Ω	±0.35 K / ±0.13 Ω	±0.27 K / ±0.10 Ω	±0.08 K / ±0.03 Ω		

Table reaction time	EPA-8, EPA-18	
t <sub>50</sub>	4.4 s	
t <sub>90</sub>	13.1 s	

### **Reaction time**

ng a

2

The mentioned times were measured by emersing a temperature sensor from room temperature into boiling water.

Electrical Connection PHARMA

#### **Electrical connection without transmitter**

#### With 1x M12 plug

3

# 4


#### Configuration 1st M12 plug



#### With M12 plug




**Electrical connection with transmitter** 



1: + supply 2: - supply 4...20 mA

3: not connected 4: not connected



## With 2x M12 plug




#### Configuration 2nd M12 plug



### With cable gland


#### Configuration strip terminal



#### Electrical connection with two transmitter (TFP-642, -842)

#### With 1x M12 plug (sensor 1 + sensor 2)

#### Configuration M12 plug



1: + supply (sensor 1)

2: - supply 4...20 mA (sensor 1) 3: - supply 4...20 mA (sensor 2)

4: + supply (sensor 2)

#### With fixed cable



## Fixed cable connection with 1x Pt100

wh ye bn gn standard rd rd wh wh PTFE

#### With 2x M12 plug (sensor 1)

#### Configuration M12 plug



1: + supply (sensor 1)


2: - supply 4...20 mA (sensor 1)

3: not connected

4: not connected

#### With 2x M12 plug (sensor 2)

#### Configuration M12 plug





1: + supply (sensor 2)


2: - supply 4...20 mA (sensor 2)


3: not connected


4: not connected


























#### **Mechanical Connection / Installation**



 The sensors are only for use with pharmaceutical buildin system PHARMadapt EPA.

#### Transport / Storage



- No outdoor storage
- · Dry and dust free
- · Not exposed to corrosive media
- Protected against solar radiation
- · Avoiding mechanical shock and vibration
- · Storage temperature -55...90 °C (-67...194 °F)
- · Relative humidity maximum 98 %

#### Reshipment



- Sensors shall be clean and must not be contaminated with dangerous media!
- Use suitable transport packaging only to avoid damage of the equipment!

#### Disposal



- Electrical devices should not be disposed of with household trash. They must be recycled in accordance with national laws and regulations.
- Take the device directly to a specialized recycling company and do not use municipal collection points.

#### Note on 3-A Sanitary Standard 74-



Information on installation according to 3-A standard is available on our website:

www.anderson-negele.com/3A74.pdf

Click on the PDF icon to download the document.

#### **Mounting Advice**



6

 The suitable insertion length depends on the pipe diameter of the measuring point. See dimension tables in product information PHARMadapt EPA.

#### **Conventional Usage**



- · Not suitable for applications in explosive areas.
- Not suitable for applications in security-relevant equipments (SIL).

#### Cleaning / Maintenance



 In case of using pressure washers, dont't point nozzle directly to electrical connections!

#### **Standards and Guidelines**



You have to comply with applicable regulations and directives.

#### Note on CE



- Applicable directives:
   Electromagnetic Compatibility Directive 2014/30/EU
- Compliance with the applicable EU directives is identified by the CE label on the product.
- The operating company is responsible for complying with the guidelines applicable to the entire installation.

#### **Accessories**

PVC-cable with M12-connection made of 1.4305 (AISI 303), IP 69 K, unshielded

M12-PVC/4 PVC-cable 4-pin, length 5 m, 10 m, 25 m

PVC-cable with M12-connection, brass nickel-plated, IP 67, shielded

M12-PVC/4G PVC-cable 4-pin, length 5 m, 10 m, 25 m

**Programming adapter** 

MPU-P 9701 Programming adapter for MPU-4, MPU-H and MPU-M

Sealing ring

DRE-5 Sealing ring for EPA-8, Ø 5 x 1.5 mm, material EPDM (FDA compliant, USP class VI)
DRE-15 Sealing ring for EPA-18, Ø 15 x 1.5 mm, material EPDM (FDA compliant, USP class VI)

Order Code PHARMA

#### Order code for version with 1x Pt100 (A: 3-A compliant) TFP-641 for PHARMadapt EPA-8, connection head ø 49 mm, non-sensitive to vibrations & TFP-661 for PHARMadapt EPA-8, connection head ø 18 mm, electrical connection via M12 plug & TFP-681 for PHARMadapt EPA-8, connection head ø 18 mm, electrical connection via 2.5 m PTFE-cable @ TFP-841 for PHARMadapt EPA-18, connection head ø 49 mm, non-sensitive to vibrations (A) for PHARMadapt EPA-18, connection head ø 18 mm, electrical connection via M12 plug (A) TFP-861 **TFP-881** for PHARMadapt EPA-18, connection head ø 18 mm, electrical connection via 2.5 m PTFE-cable & Sensor Length EL for TFP-6xx in mm 010 length 10 mm 025 length 25 mm 050 length 50 mm 100 length 100 mm Sensor Length EL for TFP-8xx in mm length 20 mm 050 length 50 mm **Accuracy Class Pt100** 1/3B 1/10B Electrical Connection for TFP-641 and TFP-841 PGcable gland M16x1.5 M12 M12 plug, standard with MPU-LCD Transmitter Х without Transmitter for TFP-641 and TFP-841 MPU-4 programmable MPU-H HART-protocol MPU-LCD with display Transmitter for TFP-661 and TFP-861 MPU-M programmable Measurement Range (only for types with transmitter; not selectable at MPU-LCD) measuring range -10...40 °C -10...40 0...50 measuring range 0...50 °C TFP Ø 18 mm Successor: TSMF 0...100 measuring range 0...100 °C no longer available! New, modular and better! 0...150 measuring range 0...150 °C 0...200 measuring range 0...200 °C All advantages at anderson-negele.com special range хх...уу TFP-641/ **A**/ **MPU-4/** 025/ M12/ 0...100

#### Note

The clamp ring is not included in scope of delivery and must be ordered separately.

SRC-05 Clamp-tension ring for EPA-8, material 1.4301 (AISI 304)
SRC-10 Clamp-tension ring for EPA-18, material 1.4301 (AISI 304)

#### **Build-In systems**



Suitable build-in systems for temperature sensors TFP-641, -642, -661, -681, -841, -842, -861, -881 you will find in product information **Process Connection PHARMadapt EPA**.

