





All advantages at anderson-negele.com

Product Information TFP-90, -94, -96, -98

**PHARMA** 

# Pharma-Temperature Sensor with Fermenter Connection

## Application/Specified usage

- · Temperature measurement in plants of pharmaceutical industry
- · Temperature measurement in fermenter connection

#### **Application examples**

- · Process monitoring
- Monitoring of CIP-/SIP-process
- · Temperature monitoring

## Hygienic design/Process connection

- · Fermenter process connection, standard length 46 mm/52 mm
- · Easy sterilizable measuring point
- · CIP-/SIP-cleaning up to 150 °C
- · Product contacting materials compliant to FDA
- Sensor completely made of stainless steel

### Features/Advantages

- · Temperature sensor electro polished  $R_a \le 0.8 \mu m$  (standard)
- · Fermenter connection made of stainless steel 316L (1.4435)
- Integrated transmitter (optional)
- $\cdot$  Inspection certificate 3.1 in scope of delivery
- · Protection class IP 69 K (with electrical connection M12 plug)

### **Options/Accessories**

- · 2 x Pt100 (not retrofittable)
- · 2 x Pt100 with two transmitters (not retrofittable)
- · Programmable transmitter MPU-4 and MPU-M with output 4...20 mA, 2-wire
- · Integrated transmitter for HART-protocol
- · Programming adapter MPU-P 9701
- · Integrated transmitter MPU-LCD with display in connecting head
- · Pt100 chip with other classes of accuracy (1/3B, 1/10B)
- · Sensor tips with diameter 3 mm and 4 mm
- · Pre-assembled connecting cable for M12 plug
- · Fixed cable in other lengths and other material available
- · Calibration certificate, DKD-laboratory certificate
- · Surfaces with  $R_a \le 0.4 \mu m$  resp. 0.6  $\mu m$  available optional

#### **Authorizations**



## Temperature sensor TFP-90



## Temperature sensor TFP-98



2

Temperature sensor				
Process connection	Fermenter sleeve	DN25 with coupling nut, G1¼"		
Insertion length	standard	25 mm		
Connection length	standard	46 mm and 52 mm		
Materials	connecting head and coupling nut fermenter connection sealing ring	stainless steel 1.4301 (AISI 304) stainless steel 1.4435 (AISI 316L), electro polished, $R_a \le 0.8$ $\mu$ m, $R_a \le 0.6$ $\mu$ m and $R_a \le 0.4$ $\mu$ m optional EPDM		
Operating pressure		10 bar maximum		
Temperature ranges	ambient sensor tip	-50+85 °C -50+250 °C		
Sensing resistor	acc. to DIN EN 60751	Pt100		
Electrical connection	cable gland cable connection fixed cable 2.5 m fixed cable 2.5 m (> 90 °C)	M16 x 1,5 M12 plug 1.4301 (AISI 304), 4-pins LIYY 4 x 0.25 mm <sup>2</sup> PTFE 4 x 0.14 mm <sup>2</sup>		
Protection class		IP 69 K (with electrical connection M12 plug)		

Transmitter MPU-4, MPU-H, MPU-M				
Temperature ranges	ambient storage	-40+85 °C -55+90 °C		
Measuring ranges	MPU-4, MPU-H, MPU-M	standard: -1040 °C, 050 / 100 / 150 / 200 °C special ranges free programable		
Accuracy	input	< ±0.25 °C		
Temperature drift	zero, span	< 0.01 % / K		
Supply	MPU-4, MPU-H, MPU-M accuracy	835 V DC 0.01 % / V (reference: 12 V DC)		
Output	signal accuracy burden	analog 420 mA < $\pm$ 0.1 % of measurement range < 600 $\Omega$ (at U <sub>B</sub> = 24 V)		
Humidity	without condensation	098 %		

Accuracy classes of temperature sensors   Tolerances for Pt100 acc. to DIN EN 60751					
Pt100	A	1/3 B	1/10 B		
0°C/100Ω	±0.15 K / ±0.06 Ω	±0.10 K / ±0.04 Ω	±0.03 K / ±0.01 Ω		
100 °C / 138.5 Ω	±0.35 K / ±0.13 Ω	±0.27 K / ±0.10 Ω	±0.08 K / ±0.03 Ω		


**Electrical Connection PHARMA** 

#### **Electrical connection without transmitter**

#### With 1 x M12 plug

3

## Configuration 1st M12 plug





## **Electrical connection with transmitter**

## With M12 plug

## Configuration M12 plug



1: + supply 2: - supply 4...20 mA 3: not connected 4: not connected

#### With 2 x M12 plug


## Configuration 2nd M12 plug






## With cable gland

## Configuration strip terminal






#### With fixed cable

## Fixed cable connection with 1 x Pt100

wh ye bn gn standard rd rd wh wh PTFE Pt100

## Fixed cable connection with 2 x Pt100 (LIYY)

bn an 1st Pt100 wh ve rd bu pk gy 2nd Pt100




## Fixed cable connection with 2 x Pt100 (PTFE)

rd rd wh 1st Pt100 ye 2nd Pt100 vt vt

## Electrical connection with two transmitter (TFP-94)

## With 1 x M12 plug (sensor 1 + sensor 2)

## Configuration M12 plug




1: + supply (sensor 1) 2: - supply 4...20 mA (sensor 1)

3: - supply 4...20 mA (sensor 2)

4: + supply (sensor 2)

## With 2 x M12 plug (sensor 1)

## Configuration M12 plug



1: + supply (sensor 1)


2: - supply 4...20 mA (sensor 1)

3: not connected

4: not connected

## With 2 x M12 plug (sensor 2)

## Configuration M12 plug





1: + supply (sensor 2)


2: - supply 4...20 mA (sensor 2)


3: not connected


4: not connected













5








## Sensor tip diameter and response time

All temperature sensors are available with smaller sensor tips, to ensure a shorter response time. The below-mentioned times were measured by immersing a temperature sensor with room temperature into boiling water.

Sensor tip Ø 4 mm







## Transport/Storage



- · No outdoor storage
- · Dry and dust free
- Not exposed to corrosive media
- · Protected against solar radiation
- · Avoiding mechanical shock and vibration
- · Storage temperature -55...+90 °C
- · Relative humidity max. 98 %

## Reshipment



- Sensors shall be clean and must not be contaminated with dangerous media!
- Use suitable transport packaging only to avoid damage of the equipment!

## Conventional usage



- · Not suitable for applications in explosive areas.
- Not suitable for applications in security-relevant equipments (SIL).

## Standards and guidelines



 You have to comply with applicable regulations and directives.

## **Cleaning/Maintenance**



 In case of using pressure washers, dont't point nozzle directly to electrical connections!

## Disposal



- Electrical devices should not be disposed of with household trash. They must be recycled in accordance with national laws and regulations.
- Take the device directly to a specialized recycling company and do not use municipal collection points.

#### **Note on CE**



- Applicable directives:
  Electromagnetic Compatibility Directive 2014/30/EU
- Compliance with the applicable EU directives is identified by the CE label on the product.
- The operating company is responsible for complying with the guidelines applicable to the entire installation.

# Temperature Transmitter MPU-LCD with Display

## Application/Specified usage

- · 4...20 mA transmitter with LCD for Pt100 temperature sensor
- · For installation in temperature sensor
- · Sensor monitoring

## **Features**

- · 4-digit display with green backlight
- · Temperature measurement in °C and °F
- · Easy range select by one button
- · Lower costs for wiring because of 2-wire technology


## Note




See product information "MPU-LCD" for details.



Order Code PHARMA



## Order code for version with 2 x Pt100 TFP-90.2 (connecting head Ø 49 mm, non-sensitive design to vibrations, no transmitter possible!) TFP-94 (like TFP-90.2, but with higher connecting head Ø 57 mm, including 2 x transmitter) TFP-98.2 (connecting head Ø 18 mm, electrical connection 2.5 m PTFE-cable; no transmitter possible!) Connection length in mm (delivery incl. sealing ring) 46 52 Sensor length in mm 025 (special length on request) ххх Diameter sensor tip in mm (on request) 4 (on request) 6



## **Accessories**

## PVC-cabel with M12-connection made of 1.4305, IP 69 K, unshielded

M12-PVC / 4-5 m PVC-cable 4-pin, length 5 m M12-PVC / 4-10 m PVC-cable 4-pin, length 10 m M12-PVC / 4-25 m PVC-cable 4-pin, length 25 m

## PVC-cable with M12-connection made of brass plated, IP 69 K, shielded

M12-PVC / 4G-5 m PVC-cable 4-pin, length 5 m M12-PVC / 4G-10 m PVC-cable 4-pin, length 10 m M12-PVC / 4G-25 m PVC-cable 4-pin, length 25 m

**Programming adapter** 

MPU-P 9701 Programming adapter for MPU-4, MPU-H and MPU-M

Sealing ring for Fermenter connection, material EPDM

DRF-20

### Programming adapter MPU-P 9701



50016 / 7.3 / 2020-10-24 / AR / EU

NEGELE MESSTECHNIK GMBH Raiffeisenweg 7 87743 Egg an der Guenz

Phone +49 (0) 83 33 . 92 04 - 0 Fax +49 (0) 83 33 . 92 04 - 49 sales@anderson-negele.com

Tech. Support: support@anderson-negele.com Phone +49 (0) 83 33 . 92 04 - 720