

Product Information TSMP

PHARMA

Temperature Sensor Mini

Application/Specified usage

- · Temperature sensor in mini housing for pharma applications
- · Temperature measuring in pipes and vessels
- · Aseptic temperature process connections without product contact for inline, precise and fast measurement. Prefabricated thermowells and build-in systems avoid opening process.
- Demounting the sensor without opening the process and without electrical disconnection avoid downtime of the equipment at calibration and maintenance.

Application examples

- · Monitoring of CIP-/SIP-process
- Safe temperature measurement in hot steam and pressurized pipes
- · Temperature monitoring in pipes or vessels

Hygienic design/Process connection

- · Hygienic process connection with CLEANadapt or PHARMadapt
- · Versions available with EHEDG approval
- · Versions available to conform to 3-A Standard 74-
- · All wetted materials are FDA-conform
- · Sensor completely made of stainless steel
- · Complete overview of process connections: see order code
- · The Anderson-Negele CLEANadapt and PHARMadapt system offers a flowoptimized, hygienic and easily sterilizable installation solution for sensors.

Features/Advantages

- · High accuracy and high ambient temperature resistance
- · Customer offset and slope adjustment
- Flex hybrid mode with digital IO-Link and analog 4...20 mA
- · Process temperature range -50...250 °C / -58...482 °F

Options/Accessories

- · 2x RTD
- · Integrated transmitter
- · Programmable transmitters TTM.I and TTM.H using IO-Link
- Different RTDs (Pt100, Pt1000) and classes of accuracy (A, AA, AAA)
- · Fast response sensor tip ø 3 mm / 0.12 in
- Spacers for high process temperature up to 250 °C / 482 °F
- · Pre-assembled connecting cable for M12 plug
- · Hardwired cable in customer length and other material available
- Programmable with any IO-Link master
- · Add-On Instructions are available at www.anderson-negele.com/aoi

Modular design





Communication



Temperature sensor TSM with Tri-Clamp



Temperature sensor TSM for PHARMadapt ESP system



2

Temperature sensor		
Process connection	CLEANadapt PHARMadapt ESP G3/8" Sensor G3/8" PHARMadapt EPA Ingold (Fermenter) Tri-Clamp Thread Plain rod	M12, G1/2" Sensor with cap nut, sensor tip Ø 3 mm Sensor with cap nut, sensor tip Ø 4 mm 8, 18 46 mm, 52 mm 1/2", 3/4", DN10, 1", 1½", 2", 2½", 3" (DIN 32676) G1/4", sensor tip Ø 3 mm (DIN ISO 228)
Tightening torque	CLEANadapt M12 CLEANadapt G1/2"	10 Nm 20 Nm
Dimensions	insertion length probe diameter sensor tip diameter	02000 mm / 078.74 in 3, 4, 6, 8, 10, 12 mm / 0.12, 0.16, 0.24, 0.31, 0.39, 0.47 in 3, 4, 6 mm / 0.12, 0.16, 0.24 in, see dimensional drawings
Materials	connecting head, spacer wetted parts sealing ring PHARMadapt EPA, Ingold (Fermenter)	stainless steel 1.4301 / AISI 304 stainless steel 1.4435 / AISI 316L stainless steel 1.4404 / AISI 316L EPDM, USP Class VI, FDA 21 CFR 177.2600
Operating pressure	CLEANadapt PHARMadapt EPA, Ingold (Fermenter)	50 bar / 725 psi maximum 10 bar / 14.5 psi maximum
Process temperature	standard range	-50250 °C / -58482 °F
Resistance Temperature Detector (RTD)	accuracy classes	Class A: ±(0.15 + 0.002 × t) °C Class AA / 1/3 DIN B: ±(0.1 + 0.0017 × t) °C Class AAA / 1/10 DIN B: ±(0.03 + 0.0005 × t) °C
Electrical connection	plug connection hardwired cable hardwired cable	M12 plug 1.4301 / AISI 304 PVC LIYY 4x 0.25 mm² / AWG 23 (perm. process temp. ≤ 90 °C) PTFE 4x 0.14 mm² / AWG 26 (perm. process temp. ≤ 250 °C)
Protection class		IP 69 K (with electrical connection M12 plug)


Transmitter TTM.I, TTM.	Н			
Temperature ranges	ambient storage	-4095 °C / -40203 °F -5590 °C / -67194 °F		
Measuring ranges		standard °C: -1040, 050 / 100 / 150 / 200 °C standard °F: 0100, 0150, 0200, 30230, 0250 °F custom ranges programable		
Accuracy	input repeatability	≤ 0.1 K (at ambient ≤ 85 °C / 185 °F) ≤ 0.05 K		
Temperature drift	typical maximum	5 mK/K (at 25 °C / 77 °F) 10 mK/K (at 25°C / 77 °F)		
Adjustments	damping offset slope	0120 s ≤ ±10 K ≤ ±25 %		
Digital output	digital resolution master cycle time power supply	IO-Link 0.01 K ≥ 51.2 ms 1830 V DC according to IO-Link		
Analog output (TTM.H only)	signal accuracy temperature drift typical temperature drift max. effect of supply voltage variations maximum load resistance power supply	420 mA, 2-wire ≤ 0.05% of upper range limit 0.0005%/K (at 25 °C / 77 °F) 0.003%/K (at 25 °C / 77 °F) < 0.001% / V (at 24 V DC) R ≤ (V DC - 12 V): 0.024 A (at 25 °C / 77 °F), see diagram 1230 V DC		

Accuracy classes of temperature sensors Tolerances for Pt100 acc. to DIN EN 60751					
Pt100	Class A	Class AA / 1/3 DIN B	Class AAA / 1/10 DIN B		
0°C/100Ω	±0.15 K / ±0.06 Ω	±0.10 K / ±0.04 Ω	±0.03 K / ±0.01 Ω		
100 °C / 138.5 Ω	±0.35 K / ±0.13 Ω	±0.27 K / ±0.10 Ω	±0.08 K / ±0.03 Ω		

Accuracy classes of temperature sensors Tolerances for Pt1000 acc. to DIN EN 60751					
Pt1000	Class A	Class AA / 1/3 DIN B	Class AAA / 1/10 DIN B		
0°C/1000Ω	±0.15 K / ±0.6 Ω	±0.10 K / ±0.4 Ω	±0.03 K / ±0.1 Ω		
100 °C / 1385.1 Ω	±0.35 K / ±1.3 Ω	±0.27 K / ±1.0 Ω	±0.08 K / ±0.3 Ω		

Electrical connection without transmitter

1x RTD with M12 plug





Electrical connection with transmitter


1x RTD with M12 plug for analog operation

1: + power supply


2: - power supply 4...20 mA


3: not connected

4: not connected



2x RTD with M12 plug





1x RTD with M12 plug for IO-Link operation

1: + power supply 24 V DC

2: not connected

3: - power supply

4: IO-Link



Hardwired cable | PVC (LIYY)

1x RTD

WH YE BN GN



WH YE BN GN 1st RTD RD BU PK GY 2nd RTD





Hardwired cable | PTFE

1x RTD


RD RD WH WH


2x RTD

RD RD WH 1st RTD VT VT YE 2nd RTD




Load resistance diagram at ambient temperature 85 °C / 185 °F






Electrical connection | Head










Hardwired cable




















Process connection









CO2 | CLEANadapt G1/2"




E08 | PHARMadapt EPA-8




E18 | PHARMadapt EPA-18




M01 | PHARMadapt ESP G3/8" cap nut, ø 3 mm




M04 | Sensor G3/8" cap nut, ø 4 mm




NO1 | Plain rod




G03 | Thread G1/4", ø 3 mm




T05 | Tri-Clamp 1/2", 1/4"



Тхх | Tri-Clamp



Ixx | Ingold



Advice



Tighten the sensor only at the lower, marked in yellow spanner flat!

Tri-Clamp size				
Туре	ø D [mm / inch]			
T10	34.0 / 1.34			
TC1	50.5 / 1.99			
TC2	64.0 / 2.52			
T25	77.5 / 3.05			
TC3	91 0 / 3 58			

Dimensions table Ingold (Fermenter) Type Ingold L1 [mm / inch] L2 [mm / inch] I46 Ingold 46 46.0 / 1.81 76.0 / 2.99 I52 Ingold 52 52.0 / 2.05 82.0 / 3.23


Sensor tip diameter and response time

All temperature sensors are available with smaller sensor tips, to ensure a shorter response time. The mentioned times were measured by emersing a temperature sensor from room temperature into boiling water.

The response times given are typical measured values and may vary due to factors such as process connection, immersion length and medium.

ø 6 mm


 $t_{50} \le 1.8 \text{ s}$ $t_{90} \le 5.2 \text{ s}$ D: 8, 10, 12 mm



ø 4 mm


 $t_{50} \le 1.2 \text{ s}$

 $t_{90} \le 3.5 \text{ s}$ D: 6, 8, 10 mm





t₉₀ ≤ 2.2 s D: 6 mm



Note on 3-A Sanitary Standard 74-



Information on installation according to 3-A standard is available on our website:

www.anderson-negele.com/3A74.pdf

Click on the PDF icon to download the document.

Note on EHEDG Hygienic Standard Type EL Class I



Information on installation according to EHEDG standard is available on our

www.anderson-negele.com/EHEDG.pdf

Click on the PDF icon to download the document.

Mechanical connection/Installation



 Use Negele CLEANadapt or PHARMadapt system for safe operation of measuring point!

Transport/Storage



- · Do not store outside
- · Store in an area that is dry and dust-free
- Do not expose to corrosive media
- · Protect against solar radiation
- · Avoid mechanical shock and vibration
- · Storage temperature -55...90 °C / -67...194 °F
- Relative humidity max. 98 %

Cleaning/Maintenance



 When using a pressure washer, do not point the nozzle directly at the electrical connections.

Reshipment



- Sensors shall be clean and free of media or heatconductive paste and must not be contaminated with dangerous media!
- Use suitable transport packaging only to avoid damage of the equipment!

Conventional usage



- · Not suitable for applications in explosive areas.
- Not suitable for applications in safety-relevant system parts (SIL).

Standards and guidelines



Compliance with the applicable regulations and directives is mandatory.

Note on CE



- Applicable directives:
 Electromagnetic Compatibility Directive 2014/30/EU
- Compliance with the applicable EU directives is identified
- by the CE label on the product.

 The operating company is responsible for complying with the guidelines applicable to the entire installation.

Disposal



- Electrical devices should not be disposed of with household trash. They must be recycled in accordance with national laws and regulations.
- Take the device directly to a specialized recycling company and do not use municipal collection points.

Accessories

PVC-cable with M12 connection, brass nickel-plated, IP69K, shielded

M12-PVC/5G-8m 5 pin, length 8 m M12-PVC/5G-15m 5 pin, length 15 m M12-PVC/5G-30m 5 pin, length 30 m

M12-EVK M12 plug screw cap made of stainless steel (1.4305 / AISI 303) with o-ring

Note on IO-Link



Information on parameters and events are available on our website:

www.anderson-negele.com/iodd

Click on the IO-Link icon to open the website.

Order Code PHARMA

Order code

7

TSMP Temperatur Sensor Mini for Pharma Applications, material wetted parts 1.4435 / AISI 316L

•	• • •	•	
Proce	ss connection (A: 3-A conform, ©: EHEDG approval)	Proces	ss connection without media contact
T05	Tri-Clamp 1/2" and 3/4" ((A) and (E) only for 3/4")	G03	Thread G1/4", sensor tip ø 3 mm,
T10	Tri-Clamp DN10		spring loaded
TC1	Tri-Clamp 1" and 1½" (A) (E)	M01	PHARMadapt ESP G3/8" with cap nut,
TC2	Tri-Clamp 2" (A) (E)		sensor tip ø 3 mm, spring loaded
T25	Tri-Clamp 2½" (A) (E)	M04	Sensor G3/8" with cap nut,
TC3	Tri-Clamp 3" (A) (E)		sensor tip ø 4 mm, spring loaded
C01	CLEANadapt M12		· · · · ·
C02	CLEANadapt G1/2"		
N01	Plain rod .		
146	Ingold 46 mm (Fermenter)		
152	Ingold 52 mm (Fermenter)		
E08	PHARMadapt EPA-8 (A)		
E18	PHARMadapt EPA-18 (A)		

Spacer extension

- X Without spacer (permanent process temperature ≤ 100 °C / 212 °F)
- Short spacer (permanent process temperature ≤ 150 °C / 305 °F)
- H Long spacer (permanent process temperature ≤ 250 °C / 482 °F)

RTD type

- **0** 1x Pt100 A, 2-wire (probe length \leq 250 mm)
- 1 1x Pt100 AA, 2-wire (probe length ≤ 150 mm)
- 2 2x Pt100 A, 2-wire (probe length $\leq 250 mm$)
- 3 $2x Pt100 AA, 2-wire (probe length \le 150 mm)$
- 4 1x Pt100 A, 4-wire (probe length ≥ 50 mm)
- 5 1x Pt100 AA, 4-wire (probe length \geq 50 mm)
- 6 1x Pt100 AAA, 4-wire (probe length ≥ 50 mm)
- 7 2x Pt100 A, (3) 4-wire (probe length \geq 50 mm, 3-wire with sensor tip Ø 3 mm)
- 8 2x Pt100 AA, (3) 4-wire (probe length $\geq 50 \text{ mm}$, 3-wire with sensor tip $\emptyset 3 \text{ mm}$)
- 9 2x Pt100 AAA, 4-wire (probe length ≥ 50 mm)
- A 1x Pt1000 A, 2-wire
- B 1x Pt1000 AA, 2-wire
- C 2x Pt1000 A, 2-wire
- D 2x Pt1000 AA, 2-wire

Variable probe	length [mm]	Probe	length f	or proce	ss conne	ction [mi
10150	In steps of 5 mm, process connection N01: min. length 30 mm	G03 36 61	M01 37 59	M04 68 148	E08 10 25	E18 20 50
160500 5501000 11002000	In steps of 10 mm In steps of 50 mm In steps of 100 mm	75 93 100 105	83 97 160	198 234 238 249	50 100	
intermediate lengths	Not for G03, M01, M04, E08, E18 (Minimum order quantity: 3 pieces)	115 120 130 140 160				

Probe diameter

- o3 mm (standard for G03, M01)
- 04 4 mm (standard for MO4)
- **06** 6 mm (not for E08)
- 08 8 mm (not for T05, C01, E08, E18)
- 10 mm (not for Txx, C01, E08, E18)
- 12 mm (not for Txx, C01, E08, E18)

Sensor tip diameter, only for probe length ≥ 50 mm

- X Without reduction (standard for G03, M01, M04)
- For probe ø 6 mm
- 4 For probe Ø 6, 8, 10 mm
- For probe Ø 8, 10, 12 mm

